

Solida

O2 FLOWMETER BS-F25 with ADAPTER



Spesification ↘

- ▷ The listed Flowmeter are accuracy tested at 10% tolerance.
- ▷ The listed Flowmeter Are calibrated at specified inlet pressure, 21°C standar atmospheric pressure.
- ▷ Calibration and color may be varied by customer requirement.
- ▷ Pressure other than indicated may affect flometer accuracy.
- ▷ Gas use appropriate gas adapter/ connector for connecting to the proper gas source .
- ▷ Teperatur transport and storage -20°C - 60°C (Ohmeda style, Silbermann style, JIS style, Chemetron style, Puritan style, Drager style, MGI style, Etc.).
- ▷ Made in Thailand.
- ▷ F8ISO 13485: 2016.

Features ↘

- ▷ Product Code : 09-1-02218.
- ▷ 0-15 LPM flow rate, the maximum flow ranges limited to 60 LPM.
- ▷ Back inlet connection.
- ▷ Body made of chromed-brass.
- ▷ Flow control knob is in front of the item.
- ▷ Level indicator is green and flow indicator ball is stainless steel.
- ▷ CE marked.



Solida [®] CE

KEMENKES RI AKL 20403918270



PT. AMSAR PRIMA MANDIRI
Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340
www.amsarmedical.com info@amsarmedical.com

Telp. (021) 75674728