

# Aarashi

## MEDICAL GAS OUTLET VACUUM

### FEATURES



- Japanese standard.
- Console type.
- The main body uses high-quality copper alloy, well-manufactured.
- Insertion and extraction 10000 times without failure.
- 100% pressure and leak tested, can be used safely and securely.
- Built-in check valve, safely cut off the gas source when the probe is pulled out.
- Swivel pipe design for easy installation.
- Color-coded for easy gas service identification.
- Working pressure 4-8bar.

MEMBER OF:



Gaseslab Indonesia  
PERKUMPULAN ORGANISASI PERUSAHAAN ALAT ALAT  
PESERTA DARI LABORATORIUM INDONESIA



**AMSTAR**  
DURABLE, RELIABLE, COMMITTED

PT AMSAR PRIMA MANDIRI

📍 Ruko Bukit Dago AA6 No. 17-18,  
Rawakalong, Gunung Sindur,  
Bogor, Jawa Barat - 16340  
✉ info@amsarmedical.com

🌐 www.amsarmedical.com  
☎ (021) 75674728  
📞 0817 7519 86