

Portable Suction Unit

Model: ASU-200



Accessories:



- Blue
- Pink
- Green

Features

- Produk dapat dipindahkan dan cara penggunaan alat sangat sederhana dan mudah untuk dioperasikan
- Ukuran produk kecil, ringan sehingga mudah untuk dibawa
- Daya hisap tinggi dengan kemampuan hisap hingga - 80 kpa
- Kekuatan baterai yang sudah terisi penuh mampu bekerja selama 60 menit
- Power suply tegangan AC 100 v - 240 v dan DC 12 v dengan menggunakan acare Mobil Charger
- Botol Suction bawaan mampu menampung 800 - 1000 ml. tergantung kebutuhan dari pengguna



CE 2460



PT. AMSAR PRIMA MANDIRI

Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340

www.amsarmedical.com info@amsarmedical.com

Telp. (021) 75674728



Technical Specifications

Model	ASU-200
Power Consumption	35VA
Max Suction Pressure	-80Kpa (\approx -600 mmHg)
Exhaust Air Volume	24 L/min
Unit Dimensions (L x W x H)	190mm x 110mm x 175mm 140mm x 340mm x 210mm
Weight	1.7 Kg/2.0kg
Sound Pressure Level	< 70 dB (A)
Supply Voltage	AC 100-240V, 50/60Hz Acare car charger for DC12V
Battery Charge Time	3hr(s)
Battery Usage Time	1hr (during fully charged)
Battery	14.8V, 2.5Ah
Container	800ml/ 1000 ml
Operation & Storage Range	Operation Temperature: +10°C ~ +40°C Storage Temperature: -15°C ~ +50°C Humidity: 10% ~ 90% Non-condensing
Standard Compliance	IEC 60601-1 IEC 60601-1-2 ISO 10079-1
IEC 60601-1 Classification	Class I Equipment Type BF Applied Part IP21
Standard Accessories	Suction container, connection tube, power adapter and carrying bag
Optional Accessories	Car charger, Handy Controller (on/off) and battery



 KEMENKES RI AKL 20903913189



PT. AMSAR PRIMA MANDIRI
Ruko Bukit Dago Blok AA6 No. 17 Rawakalong, Gunung Sindur, Bogor - Indonesia 16340
www.amsarmedical.com info@amsarmedical.com
Telp. (021) 75674728