



DURABLE, RELIABLE, COMMITTED

# Katalog Produk 2025

# AKD/TKDN

---

Layanan dan Produk sesuai dengan PERMENKES NO. 4 Tahun 2016

Proud members of:



Our products are registered by:



Recognized by:





**PT. Amsar Prima Mandiri**  
 Centro Havana Blok MA No.8 Delatinos  
 Kel. Rawabuntu Kec. Serpong Kota  
 Tangerang Selatan



Admin : 0817-751-986  
 Sales Leader : 0819-9708-2023  
 info@amsarmedical.com



DURABLE. RELIABLE. COMMITTED

# PRODUK DALAM NEGERI

## AKD/TKDN

Company Profile



### OXYGEN FLOWMETER / REGULATOR

Harga E-catalog :  
**IDR. 3.055.000**

- Gas Oksigen (Kemurnian Hingga 99.5%), Indonesia Farmakope
- Tekanan awal 0.41 MPa
- Tekanan Operasi 0.15 MPa
- 0.34 MPa (model # P-15LH and P35LH)
- Volume Botol 250 cc
- Berat 425g (tidak termasuk air)
- Katup Relief 0.04 MPa

Link E-Catalog Scan Disini:



### HIGH PRESSURE GAS REGULATOR FLOW GENTLE + G-I

Harga E-catalog :  
**IDR. 4.446.000**

- Silinder Oksigen Silinder Tipe-G (Katup QDV = W22) dengan tekanan Maksimal 19.6MPa (200kgf/cm2)
- Gas Oksigen (Kemurnian Hingga 99.5%), Indonesia Farmakope
- Tekanan Operasional 0.15MPa (fixed-type)
- Tekanan Operasi Katup Relief 0.44MPa
- Pengukur tekanan 0 - 30 MPa (use no oil)
- Berat 740g

Link E-Catalog Scan Disini:



### FLOW GENTLE + P-I

Harga E-catalog :  
**IDR. 4.000.000**

- Gas Oksigen (Kemurnian Hingga 99.5%), Indonesia Farmakope
- Tekanan awal 0.39 - 0.44 MPa
- Volume Botol 250 cc
- Berat 425g (excluding refined water)
- Katup Relief 0.04MPa
- Rentang Aliran
- Steker Adaptor
- Konfigurasi Outlet Jenis Fitting (DISS 9/16 - 18 UNF untuk Oksigen)

Link E-Catalog Scan Disini:



### SUCTION UNIT DAN PATIENT HOSE YOX GL-I (511) WITH HOOK BASE

Harga E-catalog :  
**IDR. 7.111.000**

- Set Types Standard set S-511
- Tekanan Vakum Masuk -60 ~ -80KPa
- Pengaturan Tekanan Vakum 0 ~ -60kPa
- Botol Canister 1200ml

**Paket termasuk:**

- Neo controller, tutup tabung, pemegang tabung, botol tabung plastik, selang pasien dengan konektor
- Neo controller, tutup tabung, pemegang tabung, botol tabung kaca, selang pasien dengan penyambung
- Jenis Botol Canister Palastik Kaca
- Berat keseluruhan 1070g 1670g
- Total Dimensi W131xD185xH334 (mm)
- Berat Botol Canister 230g 830g
- Dimensi Botol Canister Dimensi maksimum: φ112xH232
- Adaptor Jenis indeks pin

Link E-Catalog Scan Disini:



### KOMINDO KOMINDO SUCTION UNIT DAN PATIENT HOSE YOX-401 I

Harga E-catalog :  
**IDR. 6.061.000**

- Tekanan masuk vakum -53.3 to -80.0KPa
- Rentang kendali tekanan 0 to -80.0KPa
- Indikator vakum 0 to -100.0KPa
- Kapasitas botol 1200.mm
- Indikator vakum, pengontrol, botol plastik atau botol kaca, selang pasien (dengan konektor selang)
- Jenis botol Gelas/plastik
- Berat 1430g(dengan Botol Plastik) / 2020g(Botol Kaca)
- Dimensi (L)176 x (W)153 x (H)388 mm
- Berat botol 273g(Botol Plastik) / 860g(Botol Kaca)
- Dimensi botol (H)234mm, φ112mm (Diameter Luar)
- Adapter Pin-index / Schrader

Link E-Catalog Scan Disini:



### KOMINDO MS SINGLE OUTLET

Harga E-catalog :  
**IDR. 2.439.846**

- Model : MS atau JS ( Japan Standart Pin )
- Jenis sistem outlet implan dinding, jenis paparan luar untuk lengan penyegelan, untuk kolom langit, untuk menggantung di langit-langit untuk anestesi dan untuk konektor.
- Terdapat 4 Variasi pilihan yaitu : Oxygen, Vacuum, Air dan N2O

Link E-Catalog Scan Disini:





**PT. Amsar Prima Mandiri**  
 Centro Havana Blok MA No.8 Delatinos  
 Kel. Rawabuntu Kec. Serpong Kota  
 Tangerang Selatan



Admin : 0817-751-986  
 Sales Leader : 0819-9708-2023  
 info@amsarmedical.com




DURABLE. RELIABLE. COMMITTED

# PRODUK DALAM NEGERI AKD/TKDN

Company Profile




## BEDHEAD



**BEDHEAD OVAL**

IDR 2.710.620



**BEDHEAD OVAL PRISMA  
 ALUMINIUM ALLOY 120CM**

IDR 2.886.000

## ALARM



**AREA ALARM DIGITAL 1 GAS** IDR 10.483.650



**AREA ALARM DIGITAL 2 GAS** IDR 16.812.600



**AREA ALARM DIGITAL 3 GAS** IDR 24.266.500



**AREA ALARM DIGITAL 4 GAS** IDR 30.757.500



## ZONE BOX VALVE



**ZONE BOX VALVE 1 GAS** IDR 5.908.000



**ZONE BOX VALVE 2 GAS** IDR 8.650.900




**ZONE BOX VALVE 3 GAS** IDR 12.026.800



**ZONE BOX VALVE 4 GAS** IDR 15.402.700



## COMPRESSED AIR SYSTEM



KEMENKES RI AKD 204033211038  
 5982/TKDN/IK/VII/2024

**MEDICAL GAS  
 COMPRESSED AIR SYSTEM**

IDR 481.160.000

scan brosur disini



## VACUUM SUCTION SYSTEM



KEMENKES RI AKD 204033211029  
 5984/TKDN/IK/VII/2024

**MEDICAL GAS  
 VACUUM SUCTION SYSTEM**

IDR 547.120.700

scan brosur disini



## MANIFOLD



KEMENKES RI AKD 20403220994  
 TKDN : 11925/TKDN/IK/X/2024

**AUTOMATIC CHANGOVER  
 MANIFOLD**

IDR 150.501.000

scan brosur disini



## OUTLET



CLOUD Medical Gas Outlet Tipe JISS



CLOUD Medical Gas Outlet Tipe Australian Standard (CIG)



CLOUD Medical Gas Outlet Tipe Ohmeda

**INSTALASI GAS MEDIK OUTLET  
OXYGEN / VACUUM / AIR / N2O**

Mulai dari  
**IDR 6.000.000 / TITIK**

scan brosur disini



## REFERENSI PROYEK



**RSUD MATRAMAN**

Proyek instalasi gas medis oxygen



**RSUD TAMAN SARI**

Proyek instalasi gas medis oxygen



**RSUD SOEDARSO PONTIANAK**

Proyek instalasi gas medis



**R. Perina Lt2**


**RSUD MAMPANG PRAPATAN**

Proyek instalasi gas medis dan nursecall system



**RSUD KOTA DEPOK**

Proyek instalasi mesin medical vacuum



**RSUD LEUWILIANG**

Proyek instalasi mesin medical compressed air

## KLIEN KAMI



## PT AMSAR PRIMA MANDIRI

Centro Havana Blok MA No.8 Delatinos  
Kel. Rawabuntu Kec. Serpong Kota Tangerang Selatan

(021) 75674728

info@amsarmedical.com

www.amsarmedical.com

0819-9708-2023 (sales leader)

@amsar.primamandiri

Amsar Prima Mandiri